

Anti-Ctnnd2 Polyclonal Antibody For Wb, Elisa - Q9Uqb3

Item Number: CM0013402

Introduction

Rabbit polyclonal antibody against human CTNND2, raised with recombinant protein (aa 966-1225). Validated for WB, ELISA; cross-reacts with human and rat. Suitable for neuronal development and Wnt signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CTNND2 Rabbit pAb
Remarks/Alias	GT24; NPRAP; CTNND2
Species	Human
GeneID (Human)	1501
GeneID	1501
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 966-1225 of human CTNND2 (NP_001323.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UQB3
Protein Weight	133kDa
Shipping	Ice bag