

Anti-Ctps1 Rabbit Pab For Wb, Elisa - P17812

Item Number: CM0027751



Introduction

CTPS1 Rabbit pAb is a polyclonal antibody raised in rabbit against recombinant human CTPS1 (aa402-591). Validated for WB and ELISA, it cross-reacts with human and mouse CTPS1. Ideal for studying CTP synthesis, lymphocyte proliferation, and cancer.

[Learn More](#)

Parameter	Value
Product Name	CTPS1 Rabbit pAb
Remarks/Alias	CTPS; GATD5; IMD24; GATD5A; CTPS1
Species	Human
GeneID	1503
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 402-591 of human CTPS1 (NP_001896.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P17812
Protein Weight	67kDa
Shipping	Ice bag