

Anti-Ctrl/Ctrl1 Rabbit Monoclonal Antibody - P40313

Item Number: CM0003105

Introduction

Rabbit monoclonal antibody against human CTRL (CTRL1) for WB, IF, IHC, and ELISA. Validated cross-reactivity with human, mouse, and rat. Predicted molecular weight of 28 kDa.

[Learn More](#)

Parameter	Value
Product Name	CTRL Rabbit mAb
Remarks/Alias	CTRL1; CTRL
Species	Human
GeneID	1506
Immunogen	Synthetic peptide corresponding to amino acids 165-264 of human CTRL (P40313)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40313
Protein Weight	28kDa
Shipping	Ice bag