


Anti-Ctsc Polyclonal Antibody For Wb, If/Icc, Elisa - P53634

Item Number: CM0025322



Introduction

High-quality CTSC rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects human and mouse dipeptidyl peptidase 1, a lysosomal protease involved in immune cell activation and skin barrier function.

[Learn More](#)

Parameter	Value
Product Name	CTSC Rabbit pAb
Remarks/Alias	JP; HMS; JPD; PLS; CPPI; DPP1; DPPI; PALS; DPP-I; PDON1; CTSC
Species	Human
Gene ID (Human)	1075
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-134 of human CTSC (NP_001805.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P53634
Protein Weight	15kDa/16kDa/52kDa
Shipping	Ice bag