

Anti-Ctsh Rabbit Polyclonal Antibody For Wb, Elisa - P09668

Item Number: CM0026651



Introduction

Anti-CTSH rabbit polyclonal antibody validated for Western blot and ELISA, with cross-reactivity to human, mouse, and rat. Targets pro-cathepsin H, a lysosomal cysteine protease involved in protein degradation. Ideal for studying lysosomal function and protease activity.

[Learn More](#)

Parameter	Value
Product Name	CTSH Rabbit pAb
Remarks/Alias	ACC4; ACC5; CPSB; ACC-4; ACC-5; CTSH
Species	Human
GeneID (Human)	1512
GeneID	1512
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-140 of human CTSH (NP_004381.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09668
Protein Weight	37kDa
Shipping	Ice bag