


Anti-Ctss Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P25774

Item Number: CM0021035



Introduction

Anti-CTSS Rabbit Polyclonal Antibody validated for WB, IF/ICC, and ELISA; cross-reacts with human, mouse, and rat. Suitable for MHC class II antigen presentation and itch signaling research.

[Learn More](#)

Parameter	Specification
Product Name	CTSS Rabbit pAb
Remarks/Alias	CTSS
Species	Human
Gene ID (Human)	1520
Gene ID	1520
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 115-331 of human CTSS (NP_004070.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25774
Protein Weight	37kDa
Shipping	Ice bag