

Anti-Ctsv (Cathepsin L2) Rabbit Polyclonal Antibody For Wb And Elisa - O60911

Item Number: CM0026638



Introduction

Rabbit polyclonal antibody specific to human CTSV (cathepsin L2). Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for corneal physiology and protease research.

[Learn More](#)

Parameter	Value
Product Name	CTSV Rabbit pAb
Remarks/Alias	CTSU; CATL2; CTSL2; CTSV
Species	Human
Gene ID (Human)	1515
Gene ID	1515
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 105-334 of human CTSV (NP_001324.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60911
Protein Weight	37kDa
Shipping	Ice bag