

Anti-Ctu2 Polyclonal Antibody For Wb, Elisa - Q2Vpk5

Item Number: CM0014744

Introduction

High-quality rabbit polyclonal antibody targeting human CTU2 (Cytoplasmic tRNA 2-thiolation protein 2), suitable for Western blot and ELISA. Confirmed reactivity with human and rat samples. Predicted molecular weight 56 kDa. Ideal for research on tRNA modification and translation fidelity.

[Learn More](#)

Parameter	Value
Product Name	CTU2 Rabbit pAb
Remarks/Alias	MFRG; NCS2; UPF0432; C16orf84; CTU2
Species	Human
GeneID (Human)	348180
GeneID	348180
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 186-515 of human CTU2 (NP_001012777.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q2VPK5
Protein Weight	56kDa
Shipping	Ice bag