

Anti-Cuedc1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nwm3

Item Number: CM0023877



Introduction

High-quality rabbit polyclonal antibody targeting human CUEDC1, validated for WB and ELISA with cross-reactivity to rat. Ideal for studying CUE domain proteins and proteomics. - Q9NWM3.

[Learn More](#)

Parameter	Value
Product Name	CUEDC1 Rabbit pAb
Species	Human
GeneID (Human)	404093
GeneID	404093
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 187-386 of human CUEDC1 (NP_001258804.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9NWM3
Protein Weight	42kDa
Shipping	Ice bag