

Anti-Cul4B Polyclonal Antibody For Wb, Elisa - Q13620

Item Number: CM0017218



Introduction

Rabbit polyclonal antibody against human CUL4B, validated for WB and ELISA. Detects CUL4B in human and mouse samples. Suitable for ubiquitination and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	CUL4B Rabbit pAb
Remarks/Alias	SFM2; MRXSC; CUL-4B; MRXHF2; MRXS15; CUL4B
Species	Human
GeneID (human)	8450
GeneID	8450
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-152 of human CUL4B (NP_003579.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13620
Protein Weight	104kDa
Shipping	Ice bag