

Anti-Cul5 Polyclonal Antibody For Wb, Elisa - Q93034

Item Number: CM0029594



Introduction

CUL5 rabbit polyclonal antibody validated for WB and ELISA across human, mouse, and rat. Detects cullin-5, a scaffold protein in E3 ubiquitin ligase complexes. Ideal for ubiquitination and proteasomal degradation studies.

[Learn More](#)

Parameter	Value
Product Name	CUL5 Rabbit pAb
Remarks/Alias	CUL-5; VACM1; VACM-1; CUL5
Species	Human
GeneID (Human)	8065
GeneID	8065
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CUL5 (NP_003469.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q93034
Protein Weight	91kDa
Shipping	Ice bag