


Anti-Cullin 1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13616

Item Number: CM0023203



Introduction

High-quality rabbit polyclonal antibody targeting human Cullin 1. Validated for WB, IHC-P, ELISA with reactivity to human, mouse, rat. Ideal for ubiquitin-proteasome and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	Cullin 1 Rabbit pAb
Remarks/Alias	CUL1, MGC149834, MGC149835, cullin-1, Cullin 1
Species	Human
Gene ID	8454 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 435-593 of human Cullin 1 (NP_003583.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13616
Protein Weight	90kDa
Shipping	Ice bag