


Anti-Cullin 1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13616

Item Number: CM0020066



Introduction

Rabbit polyclonal antibody targeting human Cullin-1 (CUL1), validated for WB, IHC-P, IF/ICC, and ELISA. Ideal for studying SCF ubiquitin ligase complexes in cell cycle and signal transduction.

[Learn More](#)

Parameter	Value
Product Name	Cullin 1 Rabbit pAb
Remarks/Alias	Cullin 1; CUL1; cullin-1; MGC149834; MGC149835
Species	Human
GeneID	8454
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 435-593 of human Cullin 1 (NP_003583.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13616
Protein Weight	90kDa
Shipping	Ice bag