

Anti-Cullin 2 (Cul2) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13617

Item Number: CM0008925



Introduction

Cullin 2 (CUL2) rabbit monoclonal antibody for WB, IHC-P, and ELISA. Detects ~87/89 kDa band in human, mouse, and rat samples. Essential for ubiquitination and HIF-1 α degradation studies.

[Learn More](#)

Parameter	Value
Product Name	Cullin 2 Rabbit mAb
Remarks/Alias	CUL2; cullin-2; Cullin 2
Species	Human
Gene ID (Human)	8453
Gene ID	8453
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Cullin 2 (Q13617)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13617
Protein Weight	87 kDa/89 kDa
Shipping	Ice bag