

Anti-Cullin 2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13617

Item Number: CM0026001



Introduction

High-quality rabbit polyclonal antibody targeting human Cullin 2 (CUL2), validated for WB, IF/ICC, and ELISA applications. Suitable for studying E3 ubiquitin ligase complexes, HIF regulation, and proteasomal degradation pathways.

[Learn More](#)

Parameter	Value
Product Name	Cullin 2 Rabbit pAb
Remarks/Alias	CUL2; cullin-2; Cullin 2
Species	Human
Gene ID (Human)	8453
Gene ID	8453
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-310 of human Cullin 2 (NP_001185708.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q13617
Protein Weight	87kDa
Shipping	Ice bag