

Anti-Cullin 3 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13618

Item Number: CM0002943



Introduction

Cullin 3 Rabbit monoclonal antibody (Q13618), validated in WB, IHC-P, ELISA. Reacts with human, mouse, rat. Core scaffold of cullin-RING E3 ubiquitin ligase complexes. Ideal for ubiquitination pathway research.

[Learn More](#)

Parameter	Value
Product Name	Cullin 3 Rabbit mAb
Remarks/Alias	CUL-3; PHA2E; NEDAUS; Cullin 3
Species	Human
Gene ID (Human)	8452
Gene ID	8452
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 669-768 of human Cullin 3 (Q13618)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13618
Protein Weight	89kDa
Shipping	Ice bag