

Anti-Cullin 4A Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q13619

Item Number: CM0025738



Introduction

Rabbit polyclonal antibody against Cullin 4A, validated for WB, IF/ICC, IP, ELISA. Reacts with human, mouse, and rat CUL4A, a core component of cullin-RING E3 ubiquitin ligase complexes involved in protein ubiquitination and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	Cullin 4A Rabbit pAb
Remarks/Alias	CUL4A; cullin-4A; Cullin 4A
Species	Human
Gene ID (Human)	8451
Gene ID	8451
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-300 of human Cullin 4A (NP_001008895.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13619
Protein Weight	88kDa
Shipping	Ice bag