

Anti-Cux1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P39880 / Q13948

Item Number: CM0023522

Introduction

CUX1 Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA, reacting with human, mouse, and rat CUX1, a transcription factor key in neuronal differentiation and cell cycle control, ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	CUX1 Rabbit pAb
Remarks/Alias	CDP; CUX; p75; CASP; CDP1; COY1; Clox; GDDI; p100; p110; p200; CUTL1; GOLIM6; CDP/Cut; Cux/CDP; Nbla10317; CUX1
Species	Human
Gene ID	1523
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human CUX1 (NP_001189472.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P39880/Q13948
Protein Weight	72-77kDa/147-165kDa
Shipping	Ice bag