

Anti-Cwc22 Polyclonal Antibody For Wb And Elisa - Q9Hcg8

Item Number: CM0016784



Introduction

CWC22 Rabbit pAb from CamelBio detects pre-mRNA splicing factor CWC22 homolog in human and mouse. Suitable for WB and ELISA. Perfect tool for studying spliceosome function, EJC assembly, and nonsense-mediated mRNA decay.

[Learn More](#)

Parameter	Value
Product Name	CWC22 Rabbit pAb
Remarks/Alias	NCM, fSAPb, EIF4GL, CWC22
Species	Human
Gene ID	57703
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 639-908 of human CWC22 (NP_065994.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9HCG8
Protein Weight	105kDa
Shipping	Ice bag