

Anti-Cx3Cl1 Rabbit Monoclonal Antibody For Wb, Ip, Ihc-P, Elisa - O35188

Item Number: CM0008709

Introduction

Anti-CX3CL1 Rabbit Monoclonal Antibody validated for WB, IP, IHC-P, and ELISA. Targets mouse Fractalkine (O35188), a chemokine involved in immune response, leukocyte adhesion, and neuroinflammation. Suitable for inflammation and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	CX3CL1 Rabbit mAb
Remarks/Alias	ABCD-3; CX3C; Cxc3; D8Bwg0439e; FK; Scyd1
Species	Mouse
Gene ID (Human)	6376
Gene ID	20312
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-337 of mouse CX3CL1 (NP_033168.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	O35188
Protein Weight	42 kDa
Shipping	Ice bag