

Anti-Cx3CL1 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - O35188

Item Number: CM0030174

Introduction

High-quality rabbit polyclonal antibody against mouse CX3CL1 (Fractalkine) for WB, IF-P, IHC-P, ELISA. Validated in mouse samples. Ideal for chemokine signaling, inflammation, and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	CX3CL1 Rabbit pAb
Remarks/Alias	NTN; NTT; CX3C; CX3C3; SCYD1; ABCD-3; C3Xkine; fractalkine; neurotactin
Species	Mouse
GeneID (Human)	6376
GeneID	20312
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-337 of mouse CX3CL1 (NP_033168.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	O35188
Protein Weight	42 kDa
Shipping	Ice bag