

Anti-Cx3Cr1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P49238

Item Number: CM0025141



Introduction

CamelBio offers a rabbit polyclonal antibody against human CX3CR1 (fractalkine receptor), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for inflammation, immunology, and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	CX3CR1 Rabbit pAb
Remarks/Alias	V28; CCRL1; GPR13; CMKDR1; GPRV28; CMKBRL1; CX3CR1
Species	Human
GeneID (Human)	1524
GeneID	1524
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human CX3CR1 (NP_001328.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49238
Protein Weight	40kDa
Shipping	Ice bag