


Anti-Cxcl1 Rabbit Monoclonal Antibody For Wb, Elisa - P12850

Item Number: CM0004573



Introduction

Recombinant rabbit monoclonal antibody against mouse CXCL1 for Western blot (WB) and ELISA. Detects mouse CXCL1 with high specificity; ideal for inflammation and hematopoietic research.

[Learn More](#)

Parameter	Value
Product Name	CXCL1 Rabbit mAb
Remarks/Alias	KC; Fsp; N51; gro; Gro1; Mgsa; Scyb1; CXCL1
Species	Mouse
GeneID (Human)	2919
GeneID	14825
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-96 of mouse CXCL1 (NP_032202.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P12850
Protein Weight	10kDa
Shipping	Ice bag