

Anti-Cxcl10/Ip-10 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P02778

Item Number: CM0011148



Introduction

Rabbit monoclonal antibody specific for human CXCL10/IP-10, validated for WB, IHC-P, IF/ICC, and ELISA. Ideal for studying chemokine signaling, inflammation, and immune cell recruitment in research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	CXCL10/IP-10 Rabbit mAb
Remarks/Alias	C7; IFI10; INP10; IP-10; crg-2; mob-1; SCYB10; gIP-10; CXCL10/IP-10
Species	Human
GeneID (Human)	3627
GeneID	3627
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-98 of human CXCL10/IP-10 (NP_001556.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P02778
Protein Weight	11kDa
Shipping	Ice bag