

Anti-Cxcl10/Ip-10 Polyclonal Antibody For Wb, Elisa - P02778

Item Number: CM0013508



Introduction

Rabbit polyclonal antibody targeting human CXCL10/IP-10 (C-X-C motif chemokine 10). Validated for WB and ELISA, recognizes an 11 kDa secreted pro-inflammatory cytokine involved in immune cell chemotaxis and viral response.

[Learn More](#)

Parameter	Value
Product Name	CXCL10/IP-10 Rabbit pAb
Remarks/Alias	C7; IFI10; INP10; IP-10; crg-2; mob-1; SCYB10; gIP-10; CXCL10/IP-10
Species	Human
GeneID (Human)	3627
GeneID	3627
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-98 of human CXCL10/IP-10 (NP_001556.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P02778
Protein Weight	11 kDa
Shipping	Ice bag