

Anti-Cxcl10/Ip-10 Polyclonal Antibody For Wb, If/Icc, Elisa - P02778

Item Number: CM0020465



Introduction

CamelBio's CXCL10/IP-10 rabbit polyclonal antibody (P02778) is validated for WB, IF/ICC, and ELISA. It reacts with human and mouse, targeting the pro-inflammatory chemokine involved in immune cell migration and viral response.

[Learn More](#)

Parameter	Value
Product Name	CXCL10/IP-10 Rabbit pAb
Remarks/Alias	C7; IFI10; INP10; IP-10; crg-2; mob-1; SCYB10; gIP-10; CXCL10/IP-10
Species	Human
GeneID (human)	3627
GeneID	3627
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-98 of human CXCL10/IP-10 (NP_001556.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02778
Protein Weight	11kDa
Shipping	Ice bag