

Anti-Cxcl10/Ip-10 Rabbit Monoclonal Antibody For Wb And Elisa - P02778

Item Number: CM0006444



Introduction

High-specificity rabbit monoclonal antibody targeting human CXCL10/IP-10, validated for WB and ELISA. Essential for research into chemokine signaling, inflammatory diseases, and immune cell trafficking. SWISS: P02778.

[Learn More](#)

Parameter	Value
Product Name	CXCL10/IP-10 Rabbit mAb
Remarks/Alias	C7; IFI10; INP10; IP-10; crg-2; mob-1; SCYB10; gIP-10; CXCL10/IP-10
Species	Human
Gene ID (Human)	3627
Gene ID	3627
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-98 of human CXCL10/IP-10 (NP_001556.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P02778
Protein Weight	11kDa
Shipping	Ice bag