

Anti-Cxcl11 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O14625

Item Number: CM0029098



Introduction

Rabbit polyclonal antibody targeting human CXCL11 (IP-9/I-TAC). Validated for WB, IHC-P, ELISA, with cross-reactivity to rat. Ideal for chemokine and T-cell recruitment research.

[Learn More](#)

Parameter	Value
Product Name	CXCL11 Rabbit pAb
Remarks/Alias	IP9; H174; IP-9; b-R1; I-TAC; SCYB11; SCYB9B; CXCL11
Species	Human
GeneID (Human)	6373
GeneID	6373
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-94 of human CXCL11 (NP_005400.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	O14625
Protein Weight	10kDa
Shipping	Ice bag