


Anti-Cxcl12 Polyclonal Antibody For Wb - P48061

Item Number: CM0023559



Introduction

High-affinity rabbit polyclonal antibody against human CXCL12, validated for WB, IF-P, IHC-P, and ELISA, suitable for chemokine signaling, cancer, and immunology research.

[Learn More](#)

Parameter	Value
Product Name	CXCL12 Rabbit pAb
Remarks/Alias	IRH; PBSF; SDF1; TLSF; TPAR1; SCYB12; CXCL12
Species	Human
Gene ID (Human)	6387
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 40-93 of human CXCL12 (NP_000600.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48061
Protein Weight	11kDa
Shipping	Ice bag