


Anti-Cxcl12/Sdf-1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P48061

Item Number: CM0014845



Introduction

Rabbit polyclonal antibody targeting human CXCL12/SDF-1, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for chemotaxis, cancer, and immunology studies.

[Learn More](#)

Parameter	Value
Product Name	CXCL12 Rabbit pAb
Remarks/Alias	IRH; PBSF; SDF1; TLSF; TPAR1; SCYB12; CXCL12
Species	Human
Gene ID (Human)	6387
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 40-93 of human CXCL12 (NP_000600.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48061
Protein Weight	11kDa
Shipping	Ice bag