


Anti-Cxcl4/Pf4 Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P02776

Item Number: CM0005053



Introduction

Recombinant rabbit monoclonal antibody targeting human CXCL4/PF4, validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for chemokine signaling and platelet biology research.

[Learn More](#)

Parameter	Value
Product Name	CXCL4/PF4 Rabbit mAb
Remarks/Alias	PF-4; CXCL4; SCYB4; CXCL4/PF4
Species	Human
GeneID (Human)	5196
GeneID	5196
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-101 of human CXCL4/CXCL4/PF4 (P02776).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02776
Protein Weight	11kDa
Shipping	Ice bag