


Anti-Cxcl4/Pf4 Rabbit Polyclonal Antibody For Wb, Ihc, If, Elisa - P02776

Item Number: CM0021504



Introduction

High-quality rabbit polyclonal antibody against human CXCL4/PF4, validated for WB, IHC-P, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for chemokine, inflammation, and hematopoiesis research.

[Learn More](#)

Parameter	Value
Product Name	CXCL4/PF4 Rabbit pAb
Remarks/Alias	PF-4; CXCL4; SCYB4; CXCL4/PF4
Species	Human
GeneID (Human)	5196
GeneID	5196
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-101 of human CXCL4/PF4 (NP_002610.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02776
Protein Weight	11kDa
Shipping	Ice bag