

Anti-Cxcl6 Rabbit Polyclonal Antibody (Wb, Elisa) - P80162

Item Number: CM0022319



Introduction

CXCL6 rabbit polyclonal antibody validated for Western blot and ELISA. Reacts with human and mouse CXCL6/GCP-2, enabling chemokine research and IVD assay development in inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	CXCL6 Rabbit pAb
Remarks/Alias	GCP2; CKA-3; GCP-2; SCYB6; CXCL6
Species	Human
GeneID (Human)	6372
GeneID	6372
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 38-114 of human CXCL6 (NP_002984.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P80162
Protein Weight	12kDa
Shipping	Ice bag