

Anti-Cxcl9 Rabbit Polyclonal Antibody For Wb And Elisa - Q07325

Item Number: CM0020852

Introduction

Polyclonal antibody targeting human CXCL9 (MIG), raised in rabbit. Suitable for Western blot and ELISA. Recognizes the 14kDa chemokine involved in T-cell migration and inflammatory response.

[Learn More](#)

Parameter	Value
Product Name	CXCL9 Rabbit pAb
Remarks/Alias	CMK; MIG; Humig; SCYB9; crg-10; CXCL9
Species	Human
Gene ID	4283
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50 to the C-terminus of human CXCL9 (NP_002407.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q07325
Protein Weight	14kDa
Shipping	Ice bag