

Anti-Cxcr1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P25024

Item Number: CM0023470



Introduction

High-quality CXCR1 rabbit polyclonal antibody for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets the IL-8 receptor mediating neutrophil chemotaxis. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	CXCR1 Rabbit pAb
Remarks/Alias	C-C, CD128, CD181, CKR-1, IL8R1, IL8RA, CMKAR1, IL8RBA, CDw128a, C-C-CKR-1, CXCR1
Species	Human
GeneID (Human)	3577
GeneID	3577
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 250-350 of human CXCR1 (NP_000625.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25024
Protein Weight	40kDa
Shipping	Ice bag