


Anti-Cxcr3 Rabbit Polyclonal Antibody For Wb, Elisa - P49682

Item Number: CM0017784



Introduction

Rabbit polyclonal antibody targeting human CXCR3 (CD183), a chemokine receptor involved in immune cell trafficking and angiogenesis. Suitable for Western blot and ELISA; cross-reacts with human and rat. Ideal for immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	CXCR3 Rabbit pAb
Remarks/Alias	GPR9; MigR; CD182; CD183; Mig-R; CKR-L2; CMKAR3; IP10-R; CXCR3
Species	Human
GeneID (Human)	2833
GeneID	2833
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-53 of human CXCR3 (NP_001136269.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P49682
Protein Weight	41kDa
Shipping	Ice bag