

Anti-Cxcr4 Rabbit Monoclonal Antibody For Wb, Elisa - P61073

Item Number: CM0002616



Introduction

High-quality rabbit monoclonal antibody targeting human CXCR4 (C-X-C chemokine receptor type 4, CD184). Validated for WB and ELISA. CXCR4 is a GPCR involved in chemokine signaling, HIV entry, cell migration, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	CXCR4 Rabbit mAb
Remarks/Alias	FB22; HM89; LAP3; LCR1; NPYR; WHIM; CD184; LAP-3; LESTR; NPY3R; NPYRL; WHIMS; HSY3RR; NPYY3R; WHIMS1; D2S201E; CXCR4
Species	Human
Gene ID (Human)	7852
Gene ID	7852
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 253-352 of human CXCR4 (P61073)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P61073
Protein Weight	40kDa
Shipping	Ice bag