

Anti-Cxcr5 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P32302

Item Number: CM0007938

Introduction

High-quality rabbit monoclonal antibody targeting human CXCR5 (BLR1/CD185), validated for WB, IF-P, and ELISA applications. Cross-reacts with mouse and rat. Ideal for B-cell migration and lymphoma research.

[Learn More](#)

Parameter	Value
Product Name	CXCR5 Rabbit mAb
Remarks/Alias	BLR1; CD185; MDR15; CXCR5
Species	Human
GeneID (Human)	643
GeneID	643
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 273-372 of human CXCR5 (P32302)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32302
Protein Weight	42kDa
Shipping	Ice bag