

Anti-Cxcr5/Cd185 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P32302

Item Number: CM0007571

Introduction

Anti-CXCR5/CD185 Rabbit Monoclonal Antibody (clone unreported) targeting human C-X-C chemokine receptor type 5 (BLR1, MDR15). Ideal for Western blot, flow cytometry, and ELISA applications in B-cell and lymphoma studies.

[Learn More](#)

Parameter	Value
Product Name	CXCR5/CD185 Rabbit mAb
Remarks/Alias	BLR1; CD185; MDR15; CXCR5/CD185
Species	Human
Gene ID	643
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-47 of human CXCR5 (NP_001707.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P32302
Protein Weight	42kDa
Shipping	Ice bag