

Anti-Cxcr6 Rabbit Polyclonal Antibody For Wb And Elisa - O00574

Item Number: CM0023868

Introduction

Rabbit polyclonal antibody against human CXCR6 chemokine receptor, validated for Western blot and ELISA. Suitable for studying HIV/SIV coreceptor function and immune cell trafficking in lymphoid tissues.

[Learn More](#)

Parameter	Value
Product Name	CXCR6 Rabbit pAb
Remarks/Alias	BONZO; CD186; CDw186; STRL33; TYMSTR; CXCR6
Species	Human
Gene ID (Human)	10663
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human CXCR6 (NP_006555.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O00574
Protein Weight	39KDa
Shipping	Ice bag