


Anti-Cxxc1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9P0U4

Item Number: CM0025230



Introduction

Rabbit polyclonal antibody targeting human CXXC1 (CpG-binding protein), suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for nuclear protein detection in chromatin studies.

[Learn More](#)

Parameter	Value
Product Name	CXXC1 Rabbit pAb
Remarks/Alias	CFP1, CGBP, SPP1, PCCX1, PHF18, hCGBP, ZCGPC1, HsT2645, 2410002I16Rik, 5830420C16Rik, CXXC1
Species	Human
GeneID (Human)	30827
GeneID	30827
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 401-660 of human CXXC1 (NP_001095124.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P0U4
Protein Weight	76kDa
Shipping	Ice bag