


Anti-Cyb5A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, And Elisa - P00167

Item Number: CM0025153



Introduction

Rabbit polyclonal antibody targeting human CYB5A, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Generated against the N-terminal region of human CYB5A (NP_683725.1). Ideal for studying cytochrome b5 electron transport and membrane oxygenase function.

[Learn More](#)

Parameter	Value
Product Name	CYB5A Rabbit pAb
Remarks/Alias	CYB5; MCB5; METAG; CYB5A
Species	Human
Gene ID (Human)	1528
Gene ID	1528
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human CYB5A (NP_683725.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00167
Protein Weight	15kDa
Shipping	Ice bag