

Anti-Cyb5R1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uhq9

Item Number: CM0019929

Introduction

High-quality rabbit polyclonal antibody against human CYB5R1, validated for WB, IHC-P, and ELISA. Cross-reacts with rat. Ideal for studying NADH-cytochrome b5 reductase function and lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	CYB5R1 Rabbit pAb
Remarks/Alias	B5R1; B5R2; B5R.1; NQO3A2; humb5R2; CYB5R1
Species	Human
Gene ID	51706 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-305 of human CYB5R1 (NP_057327.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UHQ9
Protein Weight	34kDa
Shipping	Ice bag