


Anti-Cyb5R3 Monoclonal Antibody For Wb And Elisa - P00387

Item Number: CM0002524



Introduction

Rabbit monoclonal antibody against human CYB5R3 (NADH-cytochrome b5 reductase 3) for WB and ELISA. Cross-reacts with human and rat samples. Useful for studying electron transport and lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	CYB5R3 Rabbit mAb
Remarks/Alias	B5R; DIA1; CYB5R3
Species	Human
GeneID (Human)	1727
GeneID	1727
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CYB5R3 (P00387).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P00387
Protein Weight	34kDa
Shipping	Ice bag