


Anti-Cyb5R3 Polyclonal Antibody For Wb, If/Icc, Elisa - P00387

Item Number: CM0030602



Introduction

CYB5R3 Rabbit Polyclonal Antibody validated for WB, IF/ICC, and ELISA in human, mouse, and rat samples. Essential for studying cytochrome b5 reduction, lipid metabolism, and red blood cell biology.

[Learn More](#)

Parameter	Value
Product Name	CYB5R3 Rabbit pAb
Remarks/Alias	B5R; DIA1; CYB5R3
Species	Human
GeneID	1727
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 125-334 of human CYB5R3 (NP_001165131.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00387
Protein Weight	34kDa
Shipping	Ice bag