

Anti-Cyba Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P13498

Item Number: CM0016829



Introduction

High-quality anti-CYBA rabbit polyclonal antibody validated for western blot (WB), immunofluorescence (IF/ICC), and ELISA. Detects human, mouse, and rat CYBA/p22phox, a 21 kDa subunit of NADPH oxidase. Ideal for oxidative stress and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	CYBA Rabbit pAb
Remarks/Alias	CGD4; p22-PHOX; CYBA
Species	Human
Gene ID (Human)	1535
Gene ID	1535
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 126-195 of human CYBA (NP_000092.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13498
Protein Weight	21kDa
Shipping	Ice bag