


Anti-Cybrd1 Polyclonal Antibody For Wb, Elisa - Q53Tn4

Item Number: CM0021951



Introduction

Rabbit polyclonal antibody against CYBRD1, validated for WB and ELISA with human samples. Ideal for iron absorption and ferrireductase research. Includes core specifications and technical tips.

[Learn More](#)

Parameter	Value
Product Name	CYBRD1 Rabbit pAb
Remarks/Alias	DCYTB; FRRS3; CYB561A2; CYBRD1
Species	Human
Gene ID (Human)	79901
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 215-286 of human CYBRD1 (NP_079119.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q53TN4
Protein Weight	32kDa
Shipping	Ice bag