


Anti-Cyc1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P08574

Item Number: CM0015153



Introduction

Rabbit polyclonal antibody against human CYC1, a catalytic core subunit of complex III in the mitochondrial electron transport chain. Validated for Western blot, IHC-P, and ELISA with reactivity in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	CYC1 Rabbit pAb
Remarks/Alias	UQCR4; MC3DN6; CYC1
Species	Human
GeneID (Human)	1537
GeneID	1537
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 85-291 of human CYC1 (NP_001907.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08574
Protein Weight	35kDa
Shipping	Ice bag