

Anti-Cyclin A1 Polyclonal Antibody For Wb, Elisa - P78396

Item Number: CM0017501

Introduction

Rabbit polyclonal antibody targeting human Cyclin A1 for WB and ELISA applications. Cross-reacts with mouse and rat. Useful for germline meiosis and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	Cyclin A1 Rabbit pAb
Remarks/Alias	CT146; Cyclin A1
Species	Human
GeneID (Human)	8900
GeneID	8900
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human CCNA1 (NP_001104515.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P78396
Protein Weight	52kDa
Shipping	Ice bag