

Anti-Cyclin A1/A2 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P20248 / P78396

Item Number: CM0007835



Introduction

Rabbit monoclonal antibody targeting Cyclin A1/A2, validated for WB, IP, and ELISA. Reacts with human and mouse samples. Useful for cell cycle and cancer research. P20248 / P78396.

[Learn More](#)

Parameter	Value
Product Name	Cyclin A1/A2 Rabbit mAb
Remarks/Alias	CCN1; CCNA; CT146; Cyclin A1/A2
Species	Human
Gene ID (Human)	890/8900
Immunogen	Recombinant fusion protein corresponding to amino acids 250-465 of human Cyclin A1/A2 (P20248)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P20248/P78396
Protein Weight	49 kDa/52 kDa
Shipping	Ice bag